

# Unistat® 405w

## 1-litre vacuum insulated Asahi reactor

### Requirement

This case study looks at the performance of a Unistat 405w when connected to an Asahi 1-litre vacuum insulated glass reactor and compares the performance to a reactor with no insulation.

### Method

The Unistat 405w is connected to the reactor using two 1-metre insulated metal hoses. The reactor is filled with 0.75 litre of "M90.055.03", a silicon based HTF. The Unistat 405w was connected in its 3-phase option increasing the available heating power from 1.5 kW to 3 kW.

### Results

The set-point is changed from 20 °C to 180 °C. The jacket temperature rapidly ramps bringing the process temperature exactly to 180 °C in 29 minutes.

### Setup details

Unistat® 405w & 1-litre vacuum insulated Asahi reactor

- Temperature range: -45...250 °C
- Cooling power: 1.3 kW @ 250...0 °C  
0.7 kW @ -20 °C
- Heating power: 1.5 kW/3 kW
- Pump speed: 3300 rpm
- Hoses: 2x1 m; M24x1.5 (#9325)
- HTF: DW-Therm (#6479)
- Reactor: 1-litre jacketed glass reactor
- Reactor contents: 0.75 litre M90.055.03 (#6259)
- Reactor stirrer speed: 200 rpm
- Control: process



### Unistat® 405w – Asahi reactor insulated:

This graphic shows the performance of Unistat 405w working with an insulated 1-litre glass reactor. It takes 29 minutes to reach 180 °C from 20 °C.

